

SWASH® CLEANSING WIPES – FRAGRANCE FREE

GENERAL INFORMATION

| | |
|--------------------|---|
| Article number | Lo5012-48 |
| Basic UDI-DI | PP10417SUPPORT.PRODUCT57 |
| UDI-DI | 111666475986 |
| Intended use | Intended for a primary incontinence clean-up of incontinent and/or bedridden patients as an infection prevention measure. |
| Content of package | 48 wipes |
| Content of box | 12 packages |
| Manufacturer | Arion Manufacturing B.V. – The Netherlands Rutherford 11, 6422 RE Heerlen |
| Issue date | 07/2020 |
| Revision date | N.A. |

PRODUCT INFORMATION

| | |
|-------------------------------|--------------------------------------|
| Type of product | Medical device class 1 (2017/745/EC) |
| Composition cloth | Polyester & viscose |
| Production method cloth | Spunlace |
| Thickness cloth | 50 GSM |
| Box weight | 5535 gram |
| Total net waste weight cloths | ca. 145 gram |
| Total net waste weight foil | ca. 10 gram |
| pH-value | 5 - 6 |
| Dimensions cloth | 20 cm x 30 cm |
| Dimensions package | 10 cm x 26 cm |

INGREDIENTS LOTION

| | PURPOSE | % |
|--------------------------|--|--------|
| Aqua | Solvent in which the other ingredients are dissolved. | 75-100 |
| Propylene Glycol | Water-attracting ingredient that reduces dehydration of the skin. Next to that, it is a solvent and it determines the viscosity of the lotion. | 1-5 |
| Glycerin | Water-attracting ingredient that reduces dehydration of the skin. Also used as a solvent in many cosmetic products like shampoos. | 1-5 |
| Caprylyl Glycol | Skin conditioner. Also has preservative properties. | 0.1-1 |
| Polysorbate 20 | Emulsifier. Binds the lotion and stabilizes the emulsion. | 0.1-1 |
| Hydrolyzed Jojoba Esters | Water soluble version of jojoba (<i>Simmondsia chinensis</i>) oil that improves skin hydration, reduces skin redness and helps with skin barrier recovery. | |
| Glyceryl Caprylate | Humectant. Also has preservative properties. | 0.1-1 |
| Potassium Sorbate | Natural preservative. | 0.1-1 |
| Dipropylene Glycol | Solvent for the preservative system. | 0.1-1 |

| | | |
|-------------|---|----------|
| Citric Acid | Natural ingredient used to stabilize the lotion pH at a skin-neutral level. | 0.01-0.1 |
|-------------|---|----------|

| | | |
|----------------------------|---|--|
| INSTRUCTION FOR USE | | |
| | Wet wipes for mild and frequent cleansing (e.g. for incontinence) | |
| | For single use only. Do not reuse. | |
| | Only suitable for external use. | |
| | Nursing staff is advised to use up a packet of Swash [®] within 1 month after opening. | |
| | Do not flush down the toilet. | |

| | | |
|----------------------|---|--|
| TEST RESULTS | | |
| Skin irritation test | Results show that the Swash [®] Cleansing Wipes lotion does not cause any irritation of the skin. There was no sign of edema or red skin during or after the test. | |